

To Detect SARS-CoV-2 RT-LAMP Coupled with Solid-state Nanopores

The researchers discuss a highly sensitive, specific, and rapid [severe acute respiratory syndrome coronavirus 2](#) (SARS-CoV-2) detection platform by coupling the reverse transcription loop-mediated isothermal amplification (RT-LAMP) with glass nanopore sensors. Moreover, the researchers found that the RT-LAMP coupled nanopore digital counting platform could detect SARS-CoV-2 with 98% diagnostic sensitivity and 92% diagnostic specificity.



Introduction

The [coronavirus disease 2019](#) (COVID-19) is caused by the SARS-CoV-2 infection. As of November 11, 2021, over 250 million cases of SARS-CoV-2 infection have been reported worldwide, which has resulted in over 5 million deaths.

Despite mitigation strategies and vaccination rollouts, SARS-CoV-2 continues to spread and cause reinfection in vaccinated individuals as a result of emerging variants across countries. Therefore, it is imperative to develop an accessible, early, and accurate diagnostic method to facilitate robust public [health surveillance](#) and rapid testing.

Shortcomings of the current gold-standard technique for the SARS-CoV-2 testing, the [reverse transcription-polymerase chain reaction](#) (RT-PCR), include requirements for highly trained personnel, dedicated facilities, and instrumentations.

Comparatively, the RT-LAMP method is an isothermal amplification technique that has been widely studied as an alternative for [COVID-19 testing](#). This technique serves as a potential point-of-care tool because it is a rapid, sensitive, and specific technique.

Extensive designing, validation, and optimization of the technique have been completed, coupled with diverse optical detection readout methods such as changes in turbidity caused by magnesium pyrophosphate precipitate, changes in fluorescence using dyes, colorimetric pH indicators, or gel electrophoresis followed by [UV detection](#).

While each of these methods would require a minimum concentration of the [signalling reporters](#) (c_{min}) for the readout system to distinguish between SARS-CoV-2 positive and negative samples, the concentration of the c_{min} could range from millimolar (mM) to nanomolar (nM).

Depending on the average reaction rate of the assay, the required minimal concentration can be related to the time of the reaction. Thus, reducing the detection cmin is preferred towards quick time to positive.

In an effort to achieve this, the intrinsic single-molecule sensitivity of [nanopore sensors](#) is useful, as these sensors significantly reduce the cmin. In fact, this combined approach can allow the the analyte concentration to be within the picomolar (pM) range, thus reducing the time required to determine whether the test is positive or negative.


In a previous study, a LAMP-coupled nanopore sensor for [malaria](#) nucleic acid test was demonstrated. Based on label-free electronics-based nanopore sensors, this method is a new avenue for making accessible molecular testing at the point of care. In the present study, the researchers report using the RT-LAMP coupled nanopore platform for rapid detection of the SARS-CoV-2.

The Study

Using spiked saliva samples and COVID-19-positive clinical [nasopharyngeal swab](#) samples, the researchers established detection of SARS-CoV-2 using the RT-LAMP coupled nanopore platform.

Here, the nanopore sensor detects the amplicon size and measures the concentration by the digital counting method for qualitative assessment of positive/negative results. The researchers then compared the time to determine whether the sample was positive and assess the sensitivity performances of two one-pot [RT-LAMP assays](#) targeting the nucleocapsid (N) and envelope (E) genes of SARS-CoV-2.

The researchers presented the workflow of the RT-LAMP coupled nanopore method for SARS-CoV-2 detection from sample collection, preparation, to [RNA extraction](#) from the samples. The RT-LAMP amplification reaction is performed at 65°C for 15 min. During the nanopore readout in a negative control sample, amplification is negligible, whereas in a positive case, amplicons are increased significantly, resulting in a sharp increase in event rate.



The nanopore event rate has a linear relationship with the [analyte concentration](#). The researchers reported that the event rate threshold was set at 1s⁻¹ as the criteria for a positive result. Notably, this threshold is much higher than the background event rate in the negative control (<0.029 s⁻¹), such that the false-positive rate can be minimized.


Further, the researchers validated the method by adopting two [LAMP primer](#) sets targeting the N and E genes of SARS-CoV-2, respectively. They demonstrated that the N assay was quick to reach a positive result and had better copy number sensitivity.

Therefore, they chose the N primer set for our SARS-CoV-2 [RT-LAMP assay](#) in this study. The researchers optimized the RT-LAMP assay targeting the N gene to show a limit of detection (LoD) of 65 copies at the 95% confidence level.

Significantly, the test is also highly specific against other human coronavirus targets.

The researchers performed a nanopore counting analysis on the resulting amplicons, examined by the [glass nanopore sensor](#), and benchmarked the event rate with a threshold. The method was quicker because of its intrinsic single-molecule level of sensitivity.

Further, using real samples, the researchers demonstrated that the results confirmed an excellent analytical specificity of the RT-LAMP coupled nanopore sensor against the SARS-CoV-2.



Conclusion

The current study reports an RT-LAMP coupled nanopore platform for rapid detection of SARS-CoV-2. Compared to bulk conventional [optical techniques](#), this nanopore sensor enables one to identify positive/negative samples faster, with 98% diagnostic sensitivity and 92% diagnostic specificity compared to RT-PCR.

Employing this RT-LAMP coupled electronic nanopore digital counting platform technique for [SARS-CoV-2 detection](#) has significant potential to make next-generation accessible diagnostic testing that is rapid, sensitive, and specific.

Source:

<https://www.news-medical.net/news/20211111/RT-LAMP-coupled-solid-state-nanopores-for-SARS-CoV-2-detection.aspx>